

Sample Navigator

IT'S CHEAPER

Substantial data 'that are more important' can be stored.

No limits

Less work

Overviews

High-speed inputting data

IT'S RELIABLE

- Integration with scanners and (label)printers. Avoiding of manual input prevent human errors
- Overview of Workflow at a glance.
- Follows always the correct procedure (work flow)
- Validate al your important data, including the relation of the data with the maximum of freedom of free database fields, even sets of data fields

RESELLER, Visit us at:

www.samplnavigator.com/resellers

CLEARLY DOCUMENTED

All work protocols are Visio-based.

Full storages structures

Eligible settings

LESS WORK, SAVES TIME

Sample Navigator saves time, because it helps prevent human errors, such as redoing.

Select the right work protocol in no time.

Enter the date at the keen moment they become available.

Atomize as much as possible, scan any kind of label (1d, 2d,) print unique labels, get data from other systems like ZIS and EPD automatically

Use more than 60 standard work protocol steps, make combinations, and create your own work protocols. Validate, and start working.

Sample Navigator

The best solution for handling, storing, and retrieval of your samples



LOGISTIC SAMPLE MANAGEMENT



P.O. BOX 195
NL-8200 AD Lelystad
Schoepenweg 51E
NL8243 XP Lelystad
The Netherlands
+31(0)320-750265

www.SampleNavigator.com

TUBES, LABORATORY,
STORAGE, FREEZERS,
SOFTWARE, TRACKING &
TRACING, WORKPROTOCOL,
SCANNERS, BLOOD, SEEDS,
DATAMATRIX, DNA, TISSUE,
EDTA, HEPARIN, CITRATE,
STOL, FORENSIC, CELLS,
FAECES, RACKS, PALLETS,
ALFANUMERIC, HL7, ZIS, EPD,
ORDER MANAGEMENT,

WHERE ARE YOUR SAMPLES?

SAMPLE NAVIGATOR is...

CLEAR

All your work processes are fully documented and programmed in Microsoft® Visio-files.


Knowledge of each sample in combination with the right files.

Location specified information on all your tubes.

Content specified information of every storage, freezer and container.

Use your data in the best way possible.

Users only get the work protocols they need to use.



TRACK, TRACK 'N' TRACK

SampleNavigator tracks all your samples, even when you allocated them in another tube.

TRACK the newly placed tube and see where en when it is stored.

TRACK when and by whom the samples taken out of storage

TRACK globally, or track every little detail, with the use of extensive management overviews and reports.

NO RESTRICTIONS

NOT NOW – NOT EVER

Save all data that is being processed, use customizable, fully definable fields! However, there is much more...

Create fully definable tables, such as patient records, or illness. Do not just add extra data to a sample, but also combinations of sample, data to documents, freezers, etc.

Make corrections during the workflow process, and check the data.

Do not just add data, but also use them for searches, reports, and for comparing overviews.

(Complete L.I.M.S. functionality)



HIGH SPEED PROCESSING

EVEN WITH FILLED DATABASES

SampleNavigator can handle over 25.000.000 samples in one database. Even with that much samples, you can manufacture a picklist in a blink of an eye.

Going to the next step in the work protocol is also a matter of seconds.

Storing data in a database with over 1.000.000 samples with two simultaneous users, with not cost more than a few seconds to handle.

SAFETY

Closed security access

The users can only access those function of which they have authorization for.

Furthermore, users can only review the data of which they are authorized to see. These authorizations are field based. For example: specific user can only review certain aspects of a project, like the *name* of a patient.

- Relational database with integrity control.
 - Customizable connection with the database.
-

NO MORE LOOKING UP

IN OLD FASIONED MANUALS

No more old fashioned manuals; the work protocols take over here. Fully automated, fully monitored by the computer. The work protocol will prevent the skipping of steps and making mistakes this way.

This will also make it easier for new personnel to be introduced to the workflow in the field. Less errors!
